

Leada Range



Leada Rest

Features

- Provides mild to moderate contractures of the wrist, hand, fingers and thumb
- Maintains/supports the wrist/hand/fingers in a comfortable functional position
- Maintains the fingers in extension
- "Bend To Fit" - Easily adjust to correct alignment & support joints without the use of a heat gun
- High compliance due to comfort linings and lightweight nature
- Colour coded fastenings for easy application
- Includes extra cover, wash bag and removable finger separators

Indications:

- Maintains WHF in functional position
- Hypertonic contractures
- Multiple Sclerosis
- R/A
- CVA



Removable finger separators

SIZE	WIDTH OF MCP JOINT		WRIST TO FINGER TIPS		WRIST TO ELBOW CREASE		CODE LEFT	CODE RIGHT
	INCHES	CM	INCHES	CM	INCHES	CM		
Small	2½ - 3	6 - 7.5	6 - 7	15 - 17.5	< 8	< 20	LR-S-L	LR-S-R
Medium	3 - 3½	7.5 - 8.5	7 - 8	17.5 - 20	8 - 9½	20 - 24	LR-M-L	LR-M-R
Large	3½ - 4	8.5 - 10	8 - 9½	20 - 24	9½ +	24 +	LR-L-L	LR-L-R



Leada Grip

Features

- Treats moderate to severe contractures of the wrist/hand/fingers
- Maintains hand and fingers in a functional "C" position
- "Bend To Fit" – Easily adjust to correct alignment & support joints without the use of a heat gun
- High compliance due to comfort linings and lightweight nature
- Colour coded fastenings for easy application
- Two incrementally larger palmar rolls included to gradually increase extension of the MCP and PIP joints
- Includes extra cover, wash bag and removable finger separators

Indications:

- Controls wrist drop (90° range)
- Treats severe finger contractures
- Multiple Sclerosis
- R/A
- CVA

SIZE	WRIST TO ELBOW CREASE		CODE LEFT	CODE RIGHT
	INCHES	CM		
Small	< 8	< 20	LG-S-L	LG-S-R
Medium	8 - 9½	20 - 24	LG-M-L	LG-M-R
Large	9½ +	24 +	LG-L-L	LG-L-R



Leada Air Grip

Features

- Air Therapy – Ideal for pneumatic gentle stretch therapy to the MCP/PIP joints of the fingers.
- Low Load Prolonged Stretch for moderate to severe wrist, hand, finger joint stiffness, or contractures.
- "Bend to Fit" hand mouldable orthotic base can be set at the appropriate angle to support the wrist and forearm.
- Ideal to treat adaptive tissue shortening contracture(s) of the MCP / PIP joints of the fingers of 90 degrees or more and adduction of the thumb using graduated air pressure to reverse a fistled hand.
- Provides gentle air stretch to reverse tissue shortening
- Gently opens even the most severely non-fixed contracted hands
- Colour coded straps for easy application
- Easy to apply and remove

SIZE	CODE LEFT	CODE RIGHT
Universal	LAG-L	LAG-R

Indications:

- Controls wrist drop (90° range)
- Treats severe finger contractures
- Multiple Sclerosis
- R/A
- CVA



Leada Range



Leada Multi Use Boot

Features:

- Multi use boot with Flex-E-Core technology for ease of adjustment - just bend by hand to achieve a custom fit
- Adjusts from 30° dorsiflexion to 30° plantarflexion range with side adjustment straps to maintain the desired position
- Provides recumbent heel pressure relief as the heel floats inside the splint
- Heel retaining strap and foot cover adjustment strap
- Simple to use, with soft, washable liner
- Soft walking sole for indoor use
- Includes wedge for extra toe stretch

Indications:

- Plantar Fasciitis and Achilles Tendonitis
- Night resting splint for dorsiflexion / toe-walkers
- Pressure distribution for wound management
- Drop Foot and contracture management

SIZE	CODE
Adult	LMU-A
Paediatric	LMU-P



Leada Boot

Features:

- Ankle contracture boot with heat adjustable shell
- Removable transfer pad on fleece boot
- UltraSmooth version (non-Ambulatory) to reduce sheer and friction, has non-skid bottom for assisted transfers
- Dorsi Straps provided (optional) for enhanced dorsiflexion
- Secure wrap-around hook and loop fastening
- No tools required
- Anti-rotation bar and toe rest included
- Provides recumbent heel pressure relief as the heel floats inside the boot

Indications:

- Drop Foot
- Non-fixed plantarflexion contractures
- Wound care

SIZE	CALF CIRCUMFERENCE		AMBULATORY	NON-AMBULATORY
	INCHES	CM		
Regular	< 15¾	< 40	LBM-R	LB-R
Large	> 15¾	> 40	LBM-L	LB-L